

Real-time fluorescent PCR detection of fungal plant pathogens using the Smart Cycler

Reid D. Frederick and Christine L. Snyder

USDA-ARS Foreign Disease-Weed Science Research Unit, 1301 Ditto Ave.
Ft. Detrick, MD 21702

Reprinted with permission by Cepheid
www.cepheid.com
☎ 888-838-3222  Cepheid

OBJECTIVE

Evaluate real-time PCR assays for Karnal bunt and soybean rust using the Smart Cycler (Cepheid, Inc., Sunnyvale, CA).




Figure 1: Karnal bunt: a) infected wheat plant, b) bunted kernels, and c) teliospores. Soybean rust: d) infected soybean field in Zimbabwe, e) urediniospores on soybean leaf, and f) cross-section through uredium (970 x).



Figure 2: The portable Smart Cycler TD. Real-time PCR assays can be conducted at remote/field locations using the Smart Cycler. The Smart Cycler processing block has 16 independently-controlled reaction sites (I-CORE modules) allowing different cycling parameters to be run simultaneously.

INTRODUCTION

Plant diseases greatly reduce agricultural productivity in the United States, and new pathogens are continuing to emerge that create new disease problems requiring novel control measures. Accurate and timely diagnoses of plant disease are extremely important so that appropriate controls can be implemented quickly. Disease symptoms often aid with making decisions, but a definitive diagnosis requires unambiguous pathogen identification. A major mission of the USDA/ARS Foreign Disease-Weed Science Research Unit is to examine new technologies and develop novel methods for the detection and identification of exotic plant pathogens that are deemed significant threats to United States agriculture. Karnal bunt of wheat and rust on soybeans are two fungal diseases that are currently under study.

- Karnal bunt is an important disease of international trade, and in 1996 the pathogen, *Tilletia indica* Mitra, was found for the first time in the United States in areas of Arizona and California. Another closely related smut pathogen that infects ryegrass, *T. walkeri*, produces teliospores that are morphologically similar and difficult to distinguish from teliospores of *T. indica*.
- Soybean rust is a serious disease on soybeans that results in yield reductions due to premature defoliation of infected leaves. Two different *Phakopsora* species cause soybean rust. The soybean rust pathogen in Asia and Australia is *P. pachyrhizi* Sydow, while *P. meibomia*e Arthur is the tropical Latin American or New World rust species. Soybean rust is not currently found in the continental United States; however, serious yield losses due to soybean rust are likely to occur if the pathogen gained entry into the major soybean production regions of the United States.
- Previously, real-time PCR assays that are capable of distinguishing between *T. indica* and *T. walkeri* (1) and between *P. pachyrhizi* and *P. meibomia*e (2) were developed using an ABI Prism 7700 Sequence Detection System (Applied Biosystems, Inc., Foster City, CA).

MATERIALS & METHODS

Fungal Isolates

Tilletia isolates were obtained from India, Mexico, and the United States.

Phakopsora isolates were collected from Australia, Brazil, Hawaii, India, Philippines, Puerto Rico, Taiwan, and Thailand. All fungal isolates were maintained at the USDA-ARS Biological Containment Facility at Ft. Detrick under APHIS permit. *Tilletia* isolates were grown on potato dextrose agar (PDA). Rust isolates were propagated by inoculation onto the soybean cultivar “Williams” (*P. pachyrhizi* isolates) or onto red kidney bean (*P. meibomia*ie).

DNA Extraction

Mycelia was harvested from sporidial grown potato dextrose broth cultures of the *Tilletia* isolates, and the tissue was ground in 1.5-mL microcentrifuge tubes using a plastic pestle attached to a power drill. DNA was purified using a Puregene DNA isolation kit (Gentra Systems) according to the manufacturer’s directions. *Phakopsora* urediniospores were ground in 100 µL of extraction buffer (89 mM Tris-HCl (pH 8.0), 45 mM boric acid, 0.05 mM EDTA, and 1.0% (vol/vol) β-mercaptoethanol) in microcentrifuge tubes using a plastic pestle attached to a power drill. The samples were incubated at 75 °C for 15 minutes and were centrifuged for 10 minutes at 16,000 x g to pellet debris. The supernatant was transferred to new tubes and stored at -20 °C as DNA extracts.

<i>Pme</i>	GGAAAGTAAAAGTCGTAACCAAGGTTCCGTTAGGTGAACCTGCGGAAGGATCATTAAAT	56
<i>Ppa</i>	GGAAAGTAAAAGTCGTAACCAAGGTTCCGTTAGGTGAACCTGCGGAAGGATCATTAAAT	56
<i>Pme</i>	AAAAAGCTTAAAAAGAGTGCACCTTAATTTGGTCTTAAACACAAAACCTTTTATA	112
<i>Ppa</i>	AAAAAGCT----AAAGAGTGCACCTTTATTGGGCTC- AAAACTAAACCTTTTATA	105
<i>Pme</i>	TAAACCCATTTAAACTGGCTCAGCTACTTCATTTGA--AGTT-TTTGGGCAATCAC	165
<i>Ppa</i>	TAAACCCATTT-AATTGGCTCA-----TTGATTGATAAGATCTTTGGGC-AATGGT	154
<i>Pme</i>	AGCTTTGAAAAAAGTTGCAAAATCATCTATTATATCATTGTCTATTTTTATTTTATT	221
<i>Ppa</i>	AGCTTTGAAAAAAGTTGCAAAACCCACTATTA-ATCATAATCTTTTTTTTT-----	204
<i>Pme</i>	AAATAAC-CAAAGTCACAAAGAATGTTTTATAAATTTTAAAAATATATATAACTT	276
<i>Ppa</i>	---TAACTCAAAGTCAAATAGAATGTTTTATAAATTTAAATATATATATAACTT	257
<i>Pme</i>	TTAGCAATGGATCTCTAGGCTCTCATATCGATGAAGAACACAGTGAATGTGATAA	332
<i>Ppa</i>	TTAACCAATGGATCTCTAGGCTCTCATATCGATGAAGAACACAGTGAATGTGATAA	313
	Ppa1'er →	
<i>Pme</i>	TTAATGTGAATTCGAGAATTCAGTGAATCATCAAGTTTTTGAACGCACCTGTACC	388
<i>Ppa</i>	TTAATGTGAATTCGAGAATTCAGTGAATCATCAAGTTTTTGAACGCACCTGTACC	369
	FAM probe	
<i>Pme</i>	TTTTGGTATTCAAAAGGTACACCTGTTTGGAGTGCATGAATATTCTCAAC--TC	442
<i>Ppa</i>	TTTTGGTATTCAAAAGGTACACCTGTTTGGAGTGCATGAATATTCTCAACATTA	425
	← Pme2'er	
<i>Pme</i>	TTTCTTTTATTACTTAAAG--AAAAGCATGTTGGATTTTGGAGTGCCTGTG--TT	494
<i>Ppa</i>	TTTCTTTTTT----TAAAGGGAAA--TTGTTGGATTTTGGAGTGTGCTGTGCTT	475
	← Ppa2'er	
<i>Pme</i>	TTAATATAGCTCACTTTAAATAAATAAATATATAAAA--TTCTGTATATATATAT	547
<i>Ppa</i>	TTTTTGCAGCTCACTTTAAATAAATAAATATATAAAGTTTC-AGTATAT--TTT	527
<i>Pme</i>	GGTGAATAATAACAACATTTTCATCATTAAATTTATATAAAGGAAT-ATATATAGT	602
<i>Ppa</i>	GATGTAATAATAA--AATCATTTCATCAAAAAATAAATATATG--TGAGATTTATT	580
<i>Pme</i>	ATTAATATTATT--AT-TAAATTTT----AAGACCTCAAAATCAGGTG-GACTAC	650
<i>Ppa</i>	ATAACATTA--ATTGAATGTAAATTTTTTTTAAAGACCTCAAAATCAGGTGAGACTAC	634
<i>Pme</i>	CCACTGAACCTAAGCATATCAATAAGCGGAGGA	683
<i>Ppa</i>	CCACTGAACCTAAGCATATCAATAAGCGGAGGA	667

Figure 3: Real-time PCR strategy for the amplification of diagnostic PCR products using pathogen-specific oligonucleotide primers from the ITS regions. The primer pair Ppm1 (orange box)/Ppa2 (pink box) yields a PCR product of 141 bp that is specific for *P. pachyrhizi*. The combination of PCR primers Ppm1 (orange box)/Pme2 (blue box) generates a PCR product of 138 bp that is specific for *P. meibomia*ie. Both assays use the same FAM-labeled probe (purple box).

Ti TGGCGTACTGAGATGACAGAGCTCCACTCCGAAAGCTGACAGAGTTTGGCCGAAAGCGTCTGCGACACCCGA 20
Ti TGGCGTACTGAGATGACAGAGCTCCACTCCGAAAGCTGACAGAGTTTGGCCGAAAGCGTCTGCGACACCCGA 20
Ti ATCCCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 160
Ti ATCCCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 160
Ti GCGCTGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 237
Ti GCGCTGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 240
Ti CTCACCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 317
Ti CTCACCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 320
Ti CTGGCACCAGATGAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 397
Ti CTGGCACCAGATGAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 400
Ti GTTGTCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 477
Ti GTTGTCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 480
Ti CCGCAAGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 557
Ti CCGCAAGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 560
Ti CCGCTGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 637
Ti CCGCTGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 640
Ti CGCCCAATCCGCAAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 717
Ti CGCCCAATCCGCAAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 720
Ti TCTTCACTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 797
Ti TCTTCACTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 800
Ti AAAGTCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 877
Ti AAAGTCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 880
Ti AGCCGCAAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 957
Ti AGCCGCAAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 960
Ti TAGGATGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1037
Ti TAGGATGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1040
Ti GGAGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1117
Ti GGAGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1120
Ti CTGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1197
Ti CTGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1200
Ti CAGCAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1277
Ti CAGCAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1280
Ti TCGTGAATGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1357
Ti TCGTGAATGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1360
Ti CTGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1437
Ti CTGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1440
Ti GCGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1517
Ti GCGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1520
Ti CGAAGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1597
Ti CGAAGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1600
Ti CAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1677
Ti CAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1680
Ti TTGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1757
Ti TTGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1760
Ti AGCCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1837
Ti AGCCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1840
Ti CATAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1917
Ti CATAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1920
Ti CTACAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 1997
Ti CTACAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2000
Ti ACATGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2077
Ti ACATGAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2080
Ti GAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2157
Ti GAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2160
Ti CSTAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2237
Ti CSTAGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2240
Ti CCGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2297
Ti CCGCTGAGAACACAGCTGAGTATGCTCTGAGCTGAGCTGAGCCGCTGCTGATGATGCTACTCTGCTACAGCCGCTTCC 2300




Figure 5: Kinetics of real-time PCR reactions using the Smart Cycler system. *T. indica*-specific primers Tin3/Tin10 were used in the upper panel (A), while *T. walkeri*-specific primers Tin11/Tin10 were used in the lower panel (B). Both assays were conducted in combination with a homologous internal oligonucleotide probe labeled with the fluorescent dye FAM, and fluorescence was monitored after each cycle of PCR. Representative isolates of *T. indica* and *T. walkeri* are shown.




Figure 6: Kinetics of real-time PCR reactions using the Smart Cycler system. *P. pachyrhizi*-specific primers Ppm1/Ppa2 were used in the upper panel (A), while *P. meibomiaie*-specific primers Ppm1/Pme2 were used in the lower panel (B). Both assays were conducted in combination with a homologous internal oligonucleotide probe labeled with the fluorescent dye FAM, and fluorescence was monitored after each cycle of PCR. Representative isolates of *P. pachyrhizi* and *P. meibomiaie* are shown.

Figure 4: Real-time PCR strategy for the amplification of diagnostic PCR products using pathogen-specific oligonucleotide primers from the mtDNA region. The primer pair Tin3 (yellow box)/Tin10 (green box) yields a PCR product of 212 bp that is specific for *T. indica*. The combination of PCR primers Tin11 (blue box)/Tin10 (green box) generates a PCR product of 212 bp that is specific for *T. walkeri*. Both assays use the same FAM-labeled probe (pink box).

Real-time PCR Assays:

Oligonucleotide primers specific to *T. indica*, *T. walkeri*, *P. pachyrhizi* or *P. meibomiaie* were synthesized (Life Technologies/Gibco BRL) as described (1,2). The internal oligonucleotide probes (25-mers) were labeled at the 5'-end with the fluorescent reporter dye, 6-carboxy-fluorescein (FAM) and at the 3'-end with the quencher dye, 6-carboxy-tetramethyl-rhodamine (TAMRA). The real-time PCR reactions were conducted using 25 μ L total volume and analyzed for 45 cycles using a Smart Cycler TD (Cepheid, Sunnyvale, CA) following methods recommended by the manufacturer. Optimal assay conditions were: 95 $^{\circ}$ C for 5 s, 60 $^{\circ}$ C for 60 s, and 72 $^{\circ}$ C for 5 s for the *T. indica*-specific assay and 95 $^{\circ}$ C for 15 s and 60 $^{\circ}$ C for 30 s for the *P. pachyrhizi*-specific assay.




Figure 7: Determination of the amount of template DNA required for real-time PCR using the Smart Cycler for the *T. indica*-specific (A) and the *P. pachyrhizi*-specific (B) assays. Dilutions of genomic DNA of *T. indica* isolates WL1562 and *P. pachyrhizi* isolate Taiwan '72 were used.




Figure 9: *In planta* real-time PCR assays. Soybean plants were infected with *P. pachyrhizi* isolate Taiwan '72 or *P. meibomiae* isolate Puerto Rico '76. Ten nanograms of genomic DNA from infected plants was used in each assay with *P. pachyrhizi*-specific primers Ppm1/Ppa2 (A) and *P. meibomiae*-specific primers Ppm1/Pme2 Ppa2 (B).




Figure 8: Optimization of the real-time PCR cycling conditions using the Smart Cycler for the *T. indica*-specific (A) and the *P. pachyrhizi*-specific (B) assays. Plots show the duration of the denaturation and annealing cycles (2-step), and elongation cycle (3-step).

SUMMARY

- ◆ The real-time PCR assays for *T. indica* (Karnal bunt) and *T. walkeri* (ryegrass smut) were optimized for use on the Smart Cycler using DNA extracted from mycelia from germinating teliospores.
- ◆ The real-time PCR assays for the soybean rust pathogens *P. pachyrhizi* and *P. meibomiae* were optimized for use on the Smart Cycler using DNA extracted from ungerminated urediniospores and infected soybean leaf tissue.
- ◆ Real-time PCR assays were completed in approximately 1 hour using the Smart Cycler.

REFERENCES

1. Frederick, R.D., Snyder, K.E., Tooley, P.E., Berthier-Schaad, Y., Peterson, G.B., Bonde, M.R., Schaad, N.W., and Knorr, D.A. 2000. Identification and differentiation of *Tilletia indica* and *T. walkeri* using the polymerase chain reaction. *Phytopathology* **90**: 951-960.
2. Frederick, R.D., Snyder, C.L., Peterson, G.B., and Bonde, M.R. 2000. Detection and discrimination of the soybean rust pathogens *Phakopsora pachyrhizi* and *P. meibomiae* using PCR. *Phytopathology* **90**: S25.

The United States Department of Agriculture (USDA) prohibits discrimination in its programs on the basis of race, color, national origin, sex, religion, age, disability, political beliefs, and marital or familial status. Persons with disabilities who require alternate means for communication of program information (Braille, large print, audiotape, etc.) should contact the USDA Office of Communications at (202) 720-5881 (voice) or (202) 720-7802 (TDD).

To file a complaint, write the Secretary of Agriculture, United State Department of Agriculture, Washington, DC 20250, or call (202) 720-7327 (voice) or (202) 720-1127 (TDD). The USDA is an equal employment opportunity employer.