Importing Assay Definitions

Trademark and Copyright Statements

Cepheid®, the Cepheid logo, GeneXpert®, Xpert® and Xpertise™ are trademarks of Cepheid.

Microsoft® Windows®, Windows® 7, Windows® XP and Windows® 10 are trademarks of Microsoft Corporation. All other trademarks are the property of their respective owners.

Copyright © Cepheid 2020. All rights reserved.

Disclaimers

All examples (printouts, graphics, displays, screens, etc.) are for information and illustration purposes only and shall not be used for clinical or maintenance evaluations. Data shown in sample printouts and screens do not reflect actual patient names or test results. Labels depicted in the manual may appear different from actual product labels. Cepheid makes no representations or warranties about the accuracy and reliability of the information contained in this document. The information was developed to be used by persons trained and knowledgeable in GeneXpert system operation or under the direct supervision of Cepheid technical support or service representatives. Updates to this document may be issued periodically.

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List of Figures

1 List of Figures – English

Figure 1-1	GeneXpert Dx System Window	1-2
Figure 1-2	GeneXpert Define Assays Window	1-2
Figure 1-3	Import Assay Dialog Box	1-3
Figure 1-4	DVD Drive Supplied with Infinity	1-5
Figure 1-5	Plugging in the DVD Drive	1-5
Figure 1-6	Inserting the Assay Definitions CD into the DVD Drive	1-6
Figure 1-7	Xpertise Software Home Workspace	1-7
Figure 1-8	Setup Menu	1-7
Figure 1-9	Xpertise Manage Assays Workspace	1-8
Figure 1-10	Import Assay Dialog Box	1-9

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1 Importing Assay Definitions – English

This document explains how to import assay definitions into the GeneXpert® Dx Systems and the GeneXpert Infinity Systems. To perform the procedure, use the assay definition CD provided inside the Cepheid GeneXpert assay kit.

This chapter is divided into two major sections. Select the appropriate section depending upon the model of system that you are using:

- [Section 1.2, Importing Assay Definition Files on GeneXpert Dx Systems](#)
- [Section 1.3, Importing Assay Definition Files on GeneXpert Infinity-48s and Infinity-80 Systems](#)

1.1 CD Contents

The CD provided inside the kit contains:

- Importing Assay Definitions document (this document)
- Package Insert for the assay (usually in multiple languages)
- Assay Definition Files for both GeneXpert Dx and Infinity Systems

1.2 Importing Assay Definition Files on GeneXpert Dx Systems

On GeneXpert Dx Systems, you must import the assay definition files. Depending upon your preferences, you may also need to update your LIS test codes if the system is connected to an LIS (Laboratory Information System or HIS (Hospital Information System)).

Note

Screenshots shown in this section are Windows® 7 screenshots. Windows XP and Windows 10 screens will be very similar.

Note

Although you can import assay definitions, the GeneXpert Dx software does not allow you to modify the assay definitions.

To import assay definition files:

1. Insert the assay definition CD into the computer's DVD drive.
2. In the GeneXpert Dx System window, click **Define Assays**. See [Figure 1-1](#). The Define Assays window appears. See [Figure 1-2](#).
3. In the Define Assays window, click **Import**. The Import Assay dialog box appears. See [Figure 1-3](#).


Figure 1-1. GeneXpert Dx System Window


Figure 1-2. GeneXpert Define Assays Window

- Under the **Look in:** drop-down, select the DVD drive. See [Figure 1-3](#). Under the CD directory, select the GeneXpert Systems folder to view the list of assay definition files applicable to the GeneXpert Dx software.


Figure 1-3. Import Assay Dialog Box

5. Select the assay definition (.gxa/.nxa) file inside the GeneXpert Systems folder then click the **Import** button on the Import Assay dialog box. See [Figure 1-3](#).

The new assay name and version number will appear in the assay list (on the left side of the Define Assays window) and the details about the assay appears to the right of the list as shown in [Figure 1-2](#).

6. Check the assay name and version number to verify that you have imported the correct assay definitions.
7. If you need to import additional assay definition files from the same CD, repeat [Step 4](#) through [Step 6](#).

Note

For combination assays that have multiple .gxa/.nxa files, import only the assay definition files for assays that will be performed in your lab.

8. Remove the CD from the DVD drive and store the CD in a safe location in the event it is needed in the future.
9. **(Optional)** After an assay has been imported, if your system is connected to an LIS or HIS, you must update your host test code in order to download orders to the GeneXpert Dx System and/or upload test results from the GeneXpert Dx System to the LIS or HIS. For instructions on updating host test codes, see Chapter 5 of the *GeneXpert Dx System Operator Manual*.

Note

Cepheid recommends to always confirm that LIS or HIS uploaded results match GeneXpert result after any changes to the GeneXpert or host system, including, but not limited to, changes to the following:

- GeneXpert software version
 - GeneXpert assay definition files and version
 - GeneXpert host communication settings
 - Host middleware software or configuration changes
 - LIS software or configuration changes
-

1.3 Importing Assay Definition Files on GeneXpert Infinity-48s and Infinity-80 Systems

On GeneXpert Infinity-48s and Infinity-80 Systems, you must import the assay definition files. Depending upon your preferences, you may also want to set the cartridge retention settings and update your LIS test codes if the system is connected to an LIS or HIS network.

1.3.1 Connecting the DVD Drive

The usual way to import assay definitions is by using the DVD drive provided with the system to import the assay definitions from the CD. The DVD drive must be plugged into the USB connectors on the side of the monitor.

To plug in the DVD drive and load the assay definition file CD:

1. Locate the DVD drive. The drive is supplied with the system but customers decide where to store the DVD drive. It may be located on the shelf above the waste container. Contact your lab manager if you cannot locate the DVD drive. The DVD drive supplied is shown in [Figure 1-4](#).
2. Place the DVD drive on the kiosk top.

Figure 1-4. DVD Drive Supplied with Infinity


Figure 1-5. Plugging in the DVD Drive

3. On the side of the monitor, locate the two USB connectors and plug both DVD drive connectors into the two USB connectors on the side of the monitor. See [Figure 1-5](#).
4. Press the eject button on the front of the DVD drive to open the door.
5. Insert the assay definitions CD into the DVD drive and close the DVD drive door. See [Figure 1-6](#). The green light on the front of the DVD drive will flash while the drive reads the DVD.


Figure 1-6. Inserting the Assay Definitions CD into the DVD Drive

This completes the installation of the DVD drive and the assay definitions CD. Import the assay definitions following the procedure in [Section 1.3.2](#).

1.3.2 Importing Assay Definition Files

Note

Although you can import assay definitions, the Xpertise software does not allow you to modify the assay definitions.

To import assay definition files:

1. In the Xpertise software, press the **Home** icon to display the Xpertise Software Home workspace. See [Figure 1-7](#).


Figure 1-7. Xpertise Software Home Workspace

2. In the Xpertise Software Home workspace, press the **SETUP** button. The Setup menu is displayed. See Figure 1-8.


Figure 1-8. Setup Menu

- In the Setup menu, press the **MANAGE ASSAYS** button. The Manage Assays workspace is displayed showing the list of assays previously imported. See [Figure 1-9](#).


Figure 1-9. Xpertise Manage Assays Workspace

- In the Manage Assays workspace, click **IMPORT**. The Import Assay dialog box appears. See [Figure 1-10](#).
- Under the **Look in:** drop-down, select the DVD drive. Under the CD directory, select the Infinity Systems folder to view the list of assay definition files applicable to the Xpertise software. See [Figure 1-10](#).
- Select the assay definition (.gxa/.nxa) file, then click the **OPEN** button. See [Figure 1-10](#).

The new assay name and version number will appear in the assay list in the Manage Assays workspace. See [Figure 1-9](#). Check the assay name and version number to verify that you have imported the correct assay definitions.


Figure 1-10. Import Assay Dialog Box

- If you need to import additional assay definition files from the same CD, repeat [Step 1](#) through [Step 6](#).

Note

For combination assays that have multiple .gxa/.nxa files, import only the assay definition files for assays that will be performed in your lab.

- Remove the CD from the DVD drive and store the CD in a safe location in the event it is needed in the future. Disconnect the DVD drive from the USB port connectors and store the DVD drive for future use.
- (Optional)** After an assay has been imported, if you would like to configure the automation to retain cartridges for this assay for specific conditions, see Chapter 2 in your *GeneXpert Infinity System Operator Manual*.
- (Optional)** After an assay has been imported, if your system is connected to an LIS or HIS network, you must update your host test code in order to download tests to the Infinity System and/or upload test results from the Infinity System to the LIS or HIS network. For instructions on updating host test codes, see Chapter 5 of the *GeneXpert Infinity System Operator Manual*.

Note

Cepheid recommends confirming that LIS or HIS uploaded results match GeneXpert result after any changes to the GeneXpert or host system, including, but not limited to, changes to the following:

- GeneXpert software version
 - GeneXpert assay definition files and version
 - GeneXpert host communication settings
 - Host middleware software or configuration changes
 - LIS software or configuration changes
-

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