

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50	NucleoSpin Blood (50)	Page: 1/13
Printing date: 08.06.2018	Date of issue: 06.06.2018	

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1 Product identifier

REF 740951.50
 Product name NucleoSpin Blood (50)

REACH Registration number(s): see SECTION 3.1/3.2 or
 A registration number for the substance(s) does not exist because the annual tonnage does not require registration or the substance or its use is excluded from registration.

- 1 x 15 mL B3
- 1 x 12 mL B5
- 1 x 13 mL BE
- 1 x 30 mL BW
- 1 x 1.8 mL PB
- 1 x 1-75 mg Proteinase K (lyo)

1.2 Relevant identified uses of the substance or mixture and uses advised against

Relevant identified uses
 Product for analytical use.
 Exposure Scenario Classification according REACH, RIP 3.2 Codes: SU 0-2, PC 21, PROC 15, AC 0
 The exposure scenario is integrated into sections 1-16.

Uses advised against
 not described

1.3 Details of the supplier of the safety data sheet

Manufactured by:
 MACHEREY-NAGEL GmbH & Co. KG
 Neumann-Neander-Str. 6-8, 52355 Dueren, GERMANY
 Tel.: +49 2421 969 0

E-mail: sds@mn-net.com (msds@mn-net.com)

1.4 Emergency telephone number

Outside Germany (DE): Call your regional Poisons Information Service or call local Life Saving Service.
 DE: Gemeinsames Giftinformationszentrum (GGIZ) 99089 Erfurt tel. +49 361 730 730

You find our current versions of SDS (22 languages) in Internet: <http://www.mn-net.com/SDS>

SECTION 2: Hazard identification

2.0 Classification of the complete product



GHS02 GHS07 GHS08

Signal word DANGER

Hazard identification	Hazard classes/categories
H226	Flam. Liq. 3
H302	Acute Tox. 4 oral
H315	Skin Irrit. 2
H319	Eye Irrit. 2
H334	Resp. Sens. 1
H336, H335	STOT SE 3

2.1 Classification of the substance or mixture

15 mL B3

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 2/13

Printing date: 08.06.2018

Date of issue: 06.06.2018



GHS07

Signal word

WARNING

Hazard identification

Hazard classes/categories

H302
H319

Acute Tox. 4 oral
Eye Irrit. 2

12 mL B5

Signal word

Do not need labelling as hazardous
-

No hazard class

13 mL BE

Signal word

Do not need labelling as hazardous
-

No hazard class

30 mL BW



GHS02

GHS07

Signal word

WARNING

Hazard identification

Hazard classes/categories

H226
H302
H319
H336, H335

Flam. Liq. 3
Acute Tox. 4 oral
Eye Irrit. 2
STOT SE 3

1.8 mL PB

Signal word

Do not need labelling as hazardous
-

No hazard class

1-75 mg Proteinase K (Iyo)



GHS07

GHS08

Signal word

DANGER

Hazard identification

Hazard classes/categories

H315
H319
H334

Skin Irrit. 2
Eye Irrit. 2
Resp. Sens. 1

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 3/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

2.2 Label elements

According **CLP directive** inner packages must be only labelled with GHS symbol(s) and product identifier(s) (EU 1272/2008 Annex I - 1.5.1.2).

Harmful chemicals/mixtures with signal word: **WARNING** and highly flammable chemicals/mixtures must not be labelled with H and P phrases **until 125 mL** (EU 1272/2008 Annex I - 1.5.2). This labelling exemption is NOT valid for sensiblizing substances. The irritant hazard should be eliminated, because of buffer chemicals inside.

15 mL B3



GHS07

Signal word: WARNING

12 mL B5

Do not need labelling as hazardous
Signal word: -

13 mL BE

Do not need labelling as hazardous
Signal word: -

30 mL BW



GHS02



GHS07

Signal word: WARNING

1.8 mL PB

Do not need labelling as hazardous
Signal word: -

1-75 mg Proteinase K (Iyo)



GHS07



GHS08

Signal word: DANGER

H334

May cause allergy or asthma symptoms or breathing difficulties if inhaled.

P261sh, P342+311

Avoid breathing dust/vapours. If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

2.3 Other hazards

Possible hazards from physicochemical properties

In the case of pH values are less than 5 or higher than 9 then it is irritant. Flammable properties. ---

Information pertaining to particular risks to human and possible symptoms

Cause after oral intake, impairments of health when ingested in small quantities.

Kit contains small amounts of enzymes: May cause allergy or asthma symptoms or breathing difficulties if inhaled. -

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 4/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

Information pertaining to particular risks to the environment

PBT: not applicable

vPvB: not applicable

Other hazards

SECTION 3: Composition/information on ingredients

3.1 Substances or 3.2 Mixtures

15 mL B3

Chemical:	<i>guanidine hydrochloride</i>	CAS No.:	50-01-1
Classification:	H302, Acute Tox. 4 oral, H315, Skin Irrit. 2, H319, Eye Irrit. 2		
Formula:	CH ₆ ClN ₃		
Pseudonym:	guanidinium chloride		
TSCA Inventory:	listed		
REACH Reg. No.:	01-2119977063-35-0005		
EC No.:	200-002-3	Indice No.:	607-148-00-0
RTECS:	MF4300000	MFCD:	00013026
KE No.:	KE-18111		
Concentration:	36 - <50 %		
acc. CLP (GHS):	H302, Acute Tox. 4 oral, H319, Eye Irrit. 2		

Chemical:	<i>polyoxyethylene sorbitan monolaurate (Tween® 20)</i>	CAS No.:	9005-64-5
Classification:	No criteria for classification or naming of chemical not required.		
Formula:	C ₅₈ H ₁₁₄ O ₂₆		
Pseudonym:	Sorbitan, monododecanoate, poly(oxy-1,2-ethanediyl) derivs.		
TSCA Inventory:	listed		
RTECS:	TR7400000	MFCD:	00165986
KE No.:	KE-31681		
Concentration:	10 - <20 %		
acc. CLP (GHS):	The criteria for classification are not fulfilled.		

12 mL B5

Chemical:	<i>chemicals/mixture < 1%</i>	CAS No.:	-
Classification:	No criteria for classification or naming of chemical not required.		
TSCA Inventory:	all listed, <1%		
KE No.:	listed		
Concentration:	0,1 - <1 %		
acc. CLP (GHS):	The criteria for classification are not fulfilled.		

13 mL BE

Chemical:	<i>chemicals/mixture < 1%</i>	CAS No.:	-
Classification:	No criteria for classification or naming of chemical not required.		
TSCA Inventory:	all listed, <1%		
KE No.:	listed		
Concentration:	0,1 - <1 %		
acc. CLP (GHS):	The criteria for classification are not fulfilled.		

30 mL BW

Chemical:	<i>guanidine hydrochloride</i>	CAS No.:	50-01-1
Classification:	H302, Acute Tox. 4 oral, H315, Skin Irrit. 2, H319, Eye Irrit. 2		
Formula:	CH ₆ ClN ₃		
Pseudonym:	guanidinium chloride		
TSCA Inventory:	listed		
REACH Reg. No.:	01-2119977063-35-0005		
EC No.:	200-002-3	Indice No.:	607-148-00-0
RTECS:	MF4300000	MFCD:	00013026
KE No.:	KE-18111		
Concentration:	36 - <50 %		
acc. CLP (GHS):	H302, Acute Tox. 4 oral, H319, Eye Irrit. 2		



Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 5/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

Chemical:	2-propanol	CAS No.:	67-63-0
Classification:	H225, Flam. Liq. 2, H319, Eye Irrit. 2, H336, STOT SE 3		
Formula:	C ₃ H ₈ O		
Pseudonym:	isopropanol, IPA, propan-2-ol		
TSCA Inventory:	listed		
REACH Reg. No.:	01-2119457558-25-xxxx		
EC No.:	200-661-7	Index No.:	603-117-00-0
RTECS:	NT8050000	MFCD:	00011674
KE No.:	KE-29363		
Concentration:	20 - <35 %		
acc. CLP (GHS):	H226, Flam. Liq. 3, H319, Eye Irrit. 2, H336, STOT SE 3		

1.8 mL PB

Chemical:	glycerole	CAS No.:	56-81-5
Classification:	No criteria for classification or naming of chemical not required.		
Formula:	C ₃ H ₈ O ₃		
Pseudonym:	glycerin, 1,2,3-propanetriol		
TSCA Inventory:	listed (1,2,3-Propanetriol)		
REACH Reg. No.:	01-2119471987-18-xxxx		
EC No.:	200-289-5	Index No.:	n/a
RTECS:	MA8050000	MFCD:	00004722
KE No.:	KE-29297		
Concentration:	10 - <50 %		
acc. CLP (GHS):	The criteria for classification are not fulfilled.		

1-75 mg Proteinase K (Iyo)

Chemical:	proteinase K (origin: tritirachium album)	CAS No.:	39450-01-6
Classification:	H315, Skin Irrit. 2, H319, Eye Irrit. 2, H334, Resp. Sens. 1		
Formula:	Enzyme Comm. No. 3.4.21.64, origin: tritirachium album		
TSCA Inventory:	listed (CAS 102925-54-2)		
EC No.:	254-457-8	Index No.:	647-014-00-9
RTECS:	-	MFCD:	00132129
KE No.:	not listed		
Concentration:	90 - <100 %		
acc. CLP (GHS):	H315, Skin Irrit. 2, H319, Eye Irrit. 2, H334, Resp. Sens. 1		

3.3 Remarks

When not listed, mixtures are added with water [CAS No. 7732-18-5] to 100%.

List of H and P phrases: see section 16.1

SECTION 4: First aid measures

4.1 Description of first aid measures

Place insured person out of danger zone to fresh air immediately. Ensure quiet, warmth, and provide resuscitation if necessary. If necessary contact medical advice. Take to a doctor, in a raised position if there are breathing difficulties.

4.1.1 After SKIN Contact

Remove contaminated clothing. Rinse the affected skin or mucous membrane thoroughly under running water. (If possible) use soap.

4.1.2 After EYE Contact

After contact with the eyes rinse thoroughly under running water with the eyelid wide open with eye washing bottle, eye douche or running water (protect intact eye).

4.1.3 After INHALATION of vapours

After inhalation of foam or vapour fresh air should be inhaled. Keep airways free. Administer a Dexamethasone spray as soon as possible. Ensure quiet, warmth, and provide resuscitation if necessary. In the event of respiratory distress ensure that the patient inhales oxygen. Secure the breathing, heart and circulatory function. ---

4.1.4 After ORAL Intake

After oral intake lots of water should be drunk after it has been ingested. ---

4.2 Most important symptoms and effects, both acute and delayed

Chronic Effects: May cause sensitization by skin contact, also in repeated contact of small amounts. May cause allergy or asthma symptoms or breathing difficulties if inhaled. ---

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 6/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

4.3 Indication of any immediate medical attention and special treatment needed

Inform patient respectively further measures and the possibility of long-term damages. ---

SECTION 5: Firefighting measures

5.1 Extinguishing media

Fire extinguishers appropriate to the fire classification, and, if applicable, a fire blanket must be available in a prominent location in the work area. All extinguishers like FOAM, WATER SPRAY, DRY POWDER, CARBON DIOXIDE can be used.

5.2 Special hazards arising from the substance or mixture

WARNING: Flammable (GHS regulation). May form explosive vapour-air mixtures. Formation of hazardous and caustic vapour-air mixtures possible. ---

5.3 Advice for firefighters

No, for listed product. Product package burns like paper or plastic.

5.4 Additional information

SECTION 6: Accidental release measures

6.1 Personal precautions, protective equipment and emergency procedures

Do not breathe vapours. Regular staff training is necessary, indicating hazards and precautions on the basis of operating instructions. Restrictions on activity must be observed.

6.2 Environmental precautions

not necessary

6.3 Methods and material for containment and cleaning up

Bind any escaping liquid with inert absorbent. And dispose in accordance to local regulations for the disposal of hazardous chemicals. Clean any contaminated equipment and floors with plenty of water. Collect small amounts of leaked liquid and flush with water into drains.

6.4 Reference to other sections

SECTION 7: Handling and storage

7.1 Precautions for safe handling

Handling in accordance with the test instruction, that comes with the product. Use only in well-ventilated working areas.

7.2 Conditions for safe storage, including any incompatibilities

The original product package of MACHEREY-NAGEL allows a safe storage.

Storage class (VCI): 3

Water hazard class (DE): 1

7.2.1 Requirements for stock rooms and containers

Keep original product packages tightly closed during handling and storage.

7.3 Specific end use(s)

Product for analytical use.

SECTION 8: Exposure controls/personal protection

8.1 Control parameters

15 mL B3

Chemical: *guanidine hydrochloride*

CAS No.: 50-01-1

DNEL: [inh] 3.5 mg/m³

DNEL = Derived No-Effect Level (for workers)

PNEC_(fresh water): -
PNEC = Predicted No Effect Concentration

NIOSH: not listed

[TWA] Time-weighted average to a reference period of 8 hours, [STEL] Short-term exposure limit related to a 15-minute period

OSHA: not listed

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50 NucleoSpin Blood (50) Page: 7/13
 Printing date: 08.06.2018 Date of issue: 06.06.2018

Chemical: *polyoxyethylene sorbitan monolaurate (Tween® 20)* CAS No.: 9005-64-5

12 mL B5
 Chemical: *chemicals/mixture < 1%* CAS No.: -

13 mL BE
 Chemical: *chemicals/mixture < 1%* CAS No.: -

30 mL BW
 Chemical: *guanidine hydrochloride* CAS No.: 50-01-1

DNEL: [inh] 3.5 mg/m³
 DNEL = Derived No-Effect Level (for workers)

PNEC_(fresh water): -
 PNEC = Predicted No Effect Concentration

NIOSH: not listed
 [TWA] Time-weighted average to a reference period of 8 hours, [STEL] Short-term exposure limit related to a 15-minute period

OSHA: not listed

Chemical: *2-propanol* CAS No.: 67-63-0

DNEL: [inh] 500 mg/m³
 DNEL = Derived No-Effect Level (for workers)

PNEC_(fresh water): 140.9 mg/L
 PNEC = Predicted No Effect Concentration

TRGS 900 (DE): 200 ppm / 500 mg/m³
 E/e respirable

Short-term exposure factor: 2 (II), Y
 skin resorptive (H), respiratory sensitizable (Sa), skin sensitizable (Sh), teratogenic (Z) not securely excluded / (Y) certainly excluded

SUVA(CH) MAK value: 200 ppm / 500 mg/m³

TRGS 903 (DE): [Aceton B/b, U/b] 25 mg/L
 B blood, U urine, a no limitation, b end of exposition or shift

NIOSH: [TWA] 400 ppm / 980 mg/m³

NIOSH STEL: 500 ppm / 1225 mg/m³
 [TWA] Time-weighted average to a reference period of 8 hours, [STEL] Short-term exposure limit related to a 15-minute period

OSHA: [TWA] 400 ppm / 980 mg/m³

1.8 mL PB
 Chemical: *glycerole* CAS No.: 56-81-5

DNEL: [inh] 56 mg/m³
 DNEL = Derived No-Effect Level (for workers)

PNEC_(fresh water): 0.885 mg/L
 PNEC = Predicted No Effect Concentration

TRGS 900 (DE): 200 E mg/m³
 E/e respirable

Short-term exposure factor: 2 (I), Y
 skin resorptive (H), respiratory sensitizable (Sa), skin sensitizable (Sh), teratogenic (Z) not securely excluded / (Y) certainly excluded

SUVA(CH) MAK value: 50 e* mg/m³

1-75 mg Proteinase K (Iyo)
 Chemical: *proteinase K (origin: tritirachium album)* CAS No.: 39450-01-6
 SUVA(CH) MAK value: 0,00006_{15min} mg/m³

8.2 Exposure controls

Good ventilation and extraction system in the room, floor resistant to chemicals with floor drainage and washing facilities. The highest level of cleanliness must be maintained at the workplace.

8.2.1 Respiratory protection

Use for open access of these substances for example a protection filter, class A/AX. No additional recommendations.

8.2.2 Hand protection

Yes, gloves according EN 374 (permeation time >30 min - level 2), consist of PVC, natural latex, Neopren, or Nitril (f.ex. from Ansell or KCL). Use for short times chemical resistant latex gloves with code EN 374-3 level 1.

8.2.3 Eye protection

Yes, safety glasses according EN 166 with integrated side shields or wrap-around protection.

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50	NucleoSpin Blood (50)	Page: 8/13
Printing date: 08.06.2018	Date of issue: 06.06.2018	

- 8.2.4 Skin protection**
Recommended to avoid contamination with these hazards.
- 8.2.5 Personal hygiene**
Eating, drinking, smoking, taking snuff and storage of food in work areas and at outdoor workplaces is prohibited. Avoid contact with the skin, eyes and clothing. Rinse any clothing on which the substance has been spilled, and soak it in water. Wash hands thoroughly with soap and water when stopping work and before eating, and then apply protective skin cream.

SECTION 9: Physical and chemical properties

9.1 Information on basic physical and chemical properties

15 mL B3		
Appearance: liquid	Colour: slightly yellow	Odor: odorless
pH:	3.5-4	
Specific gravity:	1.15 g/cm ³	
12 mL B5		
Appearance: liquid	Colour: colourless	Odor: odorless
pH:	7-8	
Specific gravity:	1.00 g/cm ³	
13 mL BE		
Appearance: liquid	Colour: colourless	Odor: odorless
pH:	8-9	
Specific gravity:	1.0 g/cm ³	
30 mL BW		
Appearance: liquid	Colour: colourless	Odor: alcoholic
pH:	7-8	
Flash point:	25 °C	
Specific gravity:	1.06 g/cm ³	
1.8 mL PB		
Appearance: liquid	Colour: colourless	Odor: alcoholic
Specific gravity:	1.11 g/cm ³	
1-75 mg Proteinase K (Iyo)		
Appearance: solid (Iyoph.)	Colour: slightly grey	Odor: odorless

9.2 Other information

Data for the other parameters of the mixtures are not available, because no registration and no chemical safety report is required.
Relevant Properties of Substance Group

SECTION 10: Stability and reactivity

- 10.1 Reactivity**
no further data available.
- 10.2 Chemical stability**
No known instability.
- 10.3 Possibility of hazardous reactions**
Note: Can form very reactive substances with oxidizing agents. No further data available.
- 10.4 Conditions to avoid**
Not necessary. ---
- 10.5 Incompatible materials**
Avoid contact with strong acids or alkalines.
- 10.6 Hazardous decomposition products**



Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 9/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

In the original package all parts/all reagents are safety and separated stored. Decompositions are not observed during the expiration period under recommended conditions.

SECTION 11: Toxicological information

11.1 Information on toxicological effects

Following information is valid for pure substances. Quantitative data on the toxicity of this product are not available.

15 mL B3

Chemical:	<i>guanidine hydrochloride</i>	CAS No.:	50-01-1
TSCA Inventory:	listed	California Proposition 65 List:	not listed
Australia NICNAS:	not listed	Canada CEPA 1999:	DSL yes
Japan CSCL/PRTR:	not listed, Japan PDSCL:		not listed
Japan ISHL:	not listed		
South Korea TCCA:	not listed		
Korea Exist.Chem.Inventory:	KE-18111		
LD50 _{orl rat} :	475-907 mg/kg		
LC50 _{ihl rat} :	[4h] 3181-7655 µg/m ³		
LD50 _{drm rbt} :	2000 mg/kg		

Acute Effects: Cause after oral intake, impairments of health when ingested in small quantities.

Chemical:	<i>polyoxyethylene sorbitan monolaurate (Tween® 20)</i>	CAS No.:	9005-64-5
TSCA Inventory:	listed		
Korea Exist.Chem.Inventory:	KE-31681		

12 mL B5

Chemical:	<i>chemicals/mixture < 1%</i>	CAS No.:	-
TSCA Inventory:	all listed, <1%		
Korea Exist.Chem.Inventory:	listed		

13 mL BE

Chemical:	<i>chemicals/mixture < 1%</i>	CAS No.:	-
TSCA Inventory:	all listed, <1%		
Korea Exist.Chem.Inventory:	listed		

30 mL BW

Chemical:	<i>guanidine hydrochloride</i>	CAS No.:	50-01-1
TSCA Inventory:	listed	California Proposition 65 List:	not listed
Australia NICNAS:	not listed	Canada CEPA 1999:	DSL yes
Japan CSCL/PRTR:	not listed, Japan PDSCL:		not listed
Japan ISHL:	not listed		
South Korea TCCA:	not listed		
Korea Exist.Chem.Inventory:	KE-18111		
LD50 _{orl rat} :	475-907 mg/kg		
LC50 _{ihl rat} :	[4h] 3181-7655 µg/m ³		
LD50 _{drm rbt} :	2000 mg/kg		

Acute Effects: Cause after oral intake, impairments of health when ingested in small quantities.

Chemical:	<i>2-propanol</i>	CAS No.:	67-63-0
TSCA Inventory:	listed	California Proposition 65 List:	not listed
ACGIH:	1230 ppm		
Exposure Routes:	inhalation, ingestion, skin and/or eye contact		
Target Organs:	Eyes, skin, respiratory system		
Symptoms:	irritation eyes, nose, throat; drowsiness, dizziness, headache; dry cracking skin; in animals: narcosis		
Australia NICNAS:		Canada CEPA 1999:	DSL yes
Japan CSCL/PRTR:	PAC yes, Japan PDSCL: -		
Japan ISHL:	listed ≥1,0%/≥0,1%, Article 57-2 (SDS required)		
South Korea TCCA:			
Korea Exist.Chem.Inventory:	KE-29363		
LD50 _{orl rat} :	5045 mg/kg		
LC _{LowOrl hmn} :	3570 mg/kg		
LC50 _{ihl rat} :	164h g/m ³		
LD50 _{drm rbt} :	12.8 g/kg		

www.mn-net.com

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 10/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

TRGS 905 (DE): R_F C

1.8 mL PB

Chemical: *glycerole* CAS No.: 56-81-5
 TSCA Inventory: listed (1,2,3-Propanetriol)
 Exposure Routes: inhalation, skin and/or eye contact
 Target Organs: Eyes, skin, respiratory system, kidneys
 Symptoms: irritation eyes, skin, respiratory system; headache, nausea, vomiting; kidney injury
 Japan CSCL/PRTR: not listed, Japan PDSCL: not listed
 Japan ISHL: not listed
 Korea Exist.Chem.Inventory: KE-29297
 LD50_{orl rat}: 12.6 g/kg
 LD50_{drm rbt}: >18.7 g/kg

TRGS 905 (DE): R_F C

1-75 mg Proteinase K (Iyo)

Chemical: *proteinase K (origin: tritirachium album)* CAS No.: 39450-01-6
 TSCA Inventory: listed (CAS 102925-54-2)
 Japan CSCL/PRTR: not listed
 Japan ISHL: not listed
 Korea Exist.Chem.Inventory: not listed
 Acute Effects: Cause after impairments of health when ingested in small quantities.
 Chronic Effects: May cause sensitization by skin contact, also in repeated contact of small amounts. May cause allergy or asthma symptoms or breathing difficulties if inhaled.

SECTION 12: Ecological information

12.1 Toxicity

Following information is valid for pure substances.

15 mL B3

Chemical: *guanidine hydrochloride* CAS No.: 50-01-1
 PNEC(fresh water): -
 PNEC = Predicted No Effected Concentration
 LC50_{leuciscus idus/96h}: 1759 mg/L
 LC50_{fish/96h}: [4d] 690-1850; [48h] 1758-2420 mg/L
 EC50_{daphnia/48h}: 70.2 mg/L
 EC10_{pseudomonas putita/16h}: [72h] 11.8-33.5 mg/L
 Water hazard class (DE): 1 WGK No.: 0788
 Storage class (VCI): 12

Chemical: *polyoxyethylene sorbitan monolaurate (Tween® 20)* CAS No.: 9005-64-5
 Water hazard class (DE): 1
 Storage class (VCI): 10-11

12 mL B5

Chemical: *chemicals/mixture < 1%* CAS No.: -
 Water hazard class (DE): 1
 Storage class (VCI): 12-13

13 mL BE

Chemical: *chemicals/mixture < 1%* CAS No.: -
 Water hazard class (DE): 1
 Storage class (VCI): 12-13

30 mL BW

Chemical: *guanidine hydrochloride* CAS No.: 50-01-1
 PNEC(fresh water): -
 PNEC = Predicted No Effected Concentration
 LC50_{leuciscus idus/96h}: 1759 mg/L
 LC50_{fish/96h}: [4d] 690-1850; [48h] 1758-2420 mg/L
 EC50_{daphnia/48h}: 70.2 mg/L

www.mn-net.com

Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50	NucleoSpin Blood (50)	Page: 11/13
Printing date: 08.06.2018	Date of issue: 06.06.2018	

EC10_{pseudomonas putita/16h} : [72h] 11.8-33.5 mg/L
 Water hazard class (DE): 1 WGK No.: 0788
 Storage class (VCI): 12

Chemical: *2-propanol* CAS No.: 67-63-0
 PNEC_(fresh water) : 140.9 mg/L
 PNEC = Predicted No Effect Concentration
 LC50_{fish/96h} : 1400 mg/L
 EC50_{daphnia/48h} : 13.3 g/L
 IC50_{scenedesmus quadricauda/72h} : >1000 mg/L
 EC10_{pseudomonas putita/16h} : EC5: 1050 mg/L
 Water hazard class (DE): 1 WGK No.: 0135
 Dispersion coefficient_(octanol-water) : 0.05
 Storage class (VCI): 3

1.8 mL PB

Chemical: *glycerole* CAS No.: 56-81-5
 PNEC_(fresh water) : 0.885 mg/L
 PNEC = Predicted No Effect Concentration
 LC50_{fish/96h} : >5000_{24h} mg/L
 EC50_{daphnia/48h} : >10_{24h} g/L
 IC50_{scenedesmus quadricauda/72h} : IC57_d >10 g/L
 EC10_{pseudomonas putita/16h} : EC5: >10 g/L
 Water hazard class (DE): 0
 Dispersion coefficient_(octanol-water) : -1.76
 Storage class (VCI): 10

1-75 mg Proteinase K (Iyo)

Chemical: *proteinase K (origin: tritirachium album)* CAS No.: 39450-01-6
 Water hazard class (DE): 1
 Storage class (VCI): 13

12.2 Persistence and degradability

not necessary

12.3 Bioaccumulative potential

not necessary

12.4 Mobility in soil

not necessary

12.5 Results of PBT and vPvB assessment

no data available

12.6 Other adverse effects

no additional data available

SECTION 13: Disposal considerations

Please observe local regulations for collection and disposal of hazardous waste and contact waste disposal company, where you will obtain information on laboratory waste disposal (waste code number 16 05 06).

13.1 Waste treatment methods

Normally it is possible to empty small amounts (diluted!) into drains.

SECTION 14: Transport information

UN 1993 class 3 III, **Excepted Quantities** (≤30 mL/Σ≤1 L) = ADR/ IATA E1
 or

14.1 UN number: 1993 **14.2 UN proper shipping name:** Flammable liquid, n.o.s. (2-propanol mixture)
14.3 Class: 3 **14.4 Packing group:** III

Road transport

Classification code: F1
 Limited Quantity: 5 L Tunnel restriction code: E
 Excepted Quantity: E 1 Special instructions: 640E

Air transport

PAX: 355 max. weight PAX: 60 L



Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50	NucleoSpin Blood (50)	Page: 12/13
Printing date: 08.06.2018	Date of issue: 06.06.2018	

CAO: 366 max. weight CAO: 220 L
Maritime transport
 EmS: F-E, S-E Storage category: A

14.5 Environmental hazards

none, contains only small quantities of hazardous substances

14.6 Special precautions for user

not necessary

14.7 Transport in bulk according to Annex II of MARPOL and the IBC Code

not applicable

SECTION 15: Regulatory information

15.1 Safety, health and environmental regulations/legislation specific for the substance or mixture

German act governing protection from hazardous substances (Chemicals Act / Chemikaliengesetz- ChemG), revised on August 2013
 German order governing protection from hazardous substances (Ordinance on Hazardous Substances / Gefahrstoffverordnung - GefStoffV), revised on November 2010, according to Directive 98/24/EC
 TRGS 200, German engineering rules governing the classification and labelling of hazardous substances, preparations and products, updated October 2011
 MN Leaflet/User manual, also see www.mn-net.com
 Look for your country-specific regulations.

15.2 Chemical safety assessment

not necessary for these small amounts ---

SECTION 16: Other information

16.1 List of H and P phrases

16.1.1 List of relevant H phrases

H226 Flammable liquid and vapour.
 H302 Harmful if swallowed.
 H315 Causes skin irritation.
 H319 Causes serious eye irritation.
 H334 May cause allergy or asthma symptoms or breathing difficulties if inhaled.
 H336 May cause drowsiness or dizziness.

16.1.2 List of relevant P phrases

P210 Keep away from heat/sparks/open flames/hot surfaces. No smoking.
 P260D Do not breathe vapours.
 P261sh Avoid breathing dust/vapours.
 P264W Wash with water thoroughly after handling.
 P280sh Wear protective gloves/eye protection.
 P301+312 IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
 P304+340 IF INHALED: Remove person to fresh air and keep comfortable for breathing.
 P330 Rinse mouth.
 P342+311 If experiencing respiratory symptoms: Call a POISON CENTER/doctor.

16.2 Training advice

Multiple safety training of staffs about danger and protection by using hazards in working area. Additionally training and introduction of staffs for using these products.

16.3 Recommended restriction on use

Only for professional user.
 Look about employee restrictions for young people (f. ex. 94/33/EC or DE § 22 JArbSchG)!
 Look about employee restrictions for pregnant women and nursing women (f.ex. 92/85/EEC or for DE §§ 11-13 MuSchG 2017)!
 An individual package of this product or test kit has a moderate hazardous potential.

16.4 Further information

MACHEREY-NAGEL GmbH & Co. KG provides the information contained herein in good faith being up-to-date of own realizations at revision time. This document is intended only as a guide to the appropriate precautionary handling of the material by a properly trained person using this product. Individuals receiving the information must exercise their independent judgement in determining its appropriateness for a particular purpose.
 MACHEREY-NAGEL GmbH & Co. KG makes NO REPRESENTATIONS or WARRANTIES, either expressed or implied, including without limitation any warranties of merchantability, fitness for a particular purpose with respect to the information set forth herein or the product to which the information refers. Accordingly MACHEREY-NAGEL GmbH & Co. KG will not be responsible for damages resulting

www.mn-net.com



Safety Data Sheet

according to Regulations 1907/2006/EC (REACH) and 2015/830/EU

REF: 740951.50

NucleoSpin Blood (50)

Page: 13/13

Printing date: 08.06.2018

Date of issue: 06.06.2018

from use of or reliance upon this information. See terms and conditions at the end of our price lists for additional information.

16.5 Sources of key data

Regulation 453/2010/EU REACH - REQUIREMENTS FOR THE COMPILATION OF SAFETY DATA SHEETS

Regulation 487/2013/EU, 4th adaptation of CLP regulation to technical and scientific progress

TRGS 900, German engineering rules governing limits in air at work, updated September 2016

SUVA .CH, Limits in air at work 2009, revised on 01.2009

TRGS 907, German engineering rules governing listing of substances and causes of sensitizations, updated November 2011

KÜHN, BIRETT Merkblätter Gefährliche Arbeitsstoffe (Data Sheets of Hazardous Substances)

Revisions/Updates

Reason for Revision: 2016-03 Adaptation of regulation 1221/2015/EU

2017-11 Adaption of ECHA Registration dossier